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Instructions for use 2-piece ostomy system Skin wafer Combimate, Combiplast, Fleximate, Flexiplast

General information

The protective skin wafer

The skin wafers are made of hydrocolloid, an adhesive made from natural materials such as pectin, gelatine, sodium CMC and synthetic polyisobutylene. No synthetic rubbers (elastomers) or natural rubber latex are added to the skin wafers. The skin wafers have a foam top layer to make the skin wafer more flexible, making it more comfortable to wear. The stoma opening in these skin barriers can easily be stretched if it is slightly too small **(G)**. The wafers have a large moisture absorption capacity, which is important when perspiring.

The belt connection

The Fleximate/Flexiplast wafers have a belt connection. This belt connection is made to be used with the EuroTec wide belt with 2 hooks.

Intended use

The wafers are intended for use by ostomates with Colostomy, Ileostomy or Urostomy. The wafers are intended to be used in conjunction with compatible ostomy pouches to create a system for collection of stomal output. Compatible ostomy pouches are dependent on the brand of the wafer, and are marketed as Combimate, Combiplast, Fleximate, and Flexiplast pouches.

Indications

- To be used on intact, healthy skin.
- Can also be used on sensitive or irritated skin.
- Flat skin wafers can be used for stomas on skin level.
- Convex skin wafers can be used for stomas at or below skin level, or if a flat skin barrier leads to leakage.

Contraindications

There are no contraindications for the use of these wafers.

Warnings and precautions

- Should irritation occur, discontinue use and consult your healthcare professional.
 - Do not use any oil or ointment to clean the skin, the products will not attach.
 - Do not use oil or greasy products e.g. bodylotion under the products.
 - These products are for single-use only and should not be re-used.
- Re-use may lead to increased risk of infection and cross contamination. Physical properties of the device may no longer be optimal for intended use.
- Disposal of the wafers should be done through normal household waste.

Instructions

Preparation

- Ensure all supplies are available, such as a template, curved scissors, a disposal bag, pen, wet and dry wipes, a wafer and the (already closed) compatible ostomy pouch.

Cutting the skin wafer

- Cut the included cardboard measuring card to the correct shape and size **(A)**.
- Slide the measuring card over the stoma and ensure the opening fits properly and is not too tight around the stoma. No more than 1 to 2 mm of skin must be visible around the stoma **(B)**.
- If the opening fits properly around the stoma, place the measuring card with the printed side on the skin wafer's backing paper. Trace the opening onto the skin wafer **(C)**.
- Cut the drawn shape of the stoma from the back of the skin wafer **(D)**.

Application of the skin wafer

- Clean the skin and the stoma with lukewarm water and remove soap remnants. Pat the skin dry.
- Remove the backing paper from the adhesive layer **(E)**. Touch the adhesive layer as little as possible. If you use a wafer with a plaster frame, leave the release liners of the plaster frame section in place.
- If necessary, stretch the skin wafer slightly **(H)**, to ensure the opening fits around the stoma properly.
- Apply the skin wafer to smooth skin and rub for 30 seconds, particularly around the stoma **(H)**. The skin wafer adheres best to slightly damp skin. Pre-warming is not required. The skin wafer will need some time to properly adhere to the skin. The skin wafer should be worn for no longer than 72 hours before removal.
- If you use a wafer with a plaster frame, remove the backing papers one by one and apply the plaster to the skin.

Application of the compatible pouch

All wafers have a press-to-seal system to attach the compatible pouches. Pouches with a LockRing closure (Fleximate/Flexiplast) can be additionally secured by bringing the Lockring tabs together. Wafers with a LockRing closure (Fleximate/Flexiplast) are provided with two belt tabs for additional security if a belt is required.

Removal of the wafer

Hold the skin and 'peel' the top of the skin barrier downward and off the skin in one slow motion **(I)**. Any thin layer of adhesive that remains on the skin can be covered with the next skin wafer.

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