

**INSTRUCTIONS FOR USE****HYGIENE PULP™- DISPOSABLE MID-STREAM URINE SAMPLER****REF & Description**

REF	Description
08040	Disposable Mid-Stream Urine Sampler, 1L

Product description

Disposable Mid-Stream Urine Sampler made from pulp material with water resistance capacity. Non-sterile, the product is also known as Commode Specimen Collector.

Intended use

Disposable Mid-Stream Urine Sampler is used for collecting fresh urine sample (specimen) for the purpose of analysis. Single-use, with a minimum of 4 hours water resistance.

The products are capable of retaining water and other human waste for a minimum period of 4 hours in accordance to PAS29:1999 fluid resistance test. The products have proven its capability to withstand following liquids for a maximum of 4 hours:

- Water up to 80°C
- Urine
- Sodium chloride 0.9%
- Chlorohexidine gluconate 0.05%
- Dilute Hydrochloric acid 0.1%
- Acetic acid 5%
- Hydrogen Peroxide 0.5%

The products have proven its capability to withstand following liquids, but only up to 3 hours:

- Sodium Hypochlorite 10%
- Acid (pH2-3)

Shelf Life

8 years from date of manufacture

Caution

- Do not use products if there are signs of tear.
- The products are intended for single use only. Do not reuse.
- Keep products away from direct sunlight and ultraviolet light to avoid discolouration.
- Keep products at dry place with humidity level between 30-60% RH. Avoid contact with rain, water and liquid.
- Graduated marking on product are approximate volume only.
- Do not use beyond the recommended 4 hours duration.
- Do not use on liquid that is highly acidic (pH1 to pH 3).
- Do not use for holding fat, oil, alcohol, etc..

Instruction for Use

The use and disposal of Hygiene Pulp™ should be in accordance to current acceptable clinical practice.

1. Remove the Disposable Mid-Stream Urine Sampler from the dispensing box.
2. Place it into the toilet bowl or commode chair to collect urine sample.
3. Patient will sit on the toilet bowl or commode chair to pee.
4. Remove the Urine Sampler.
5. Pour the collected urine sample into specimen cup for further analysis.
6. The contaminated Urine Sampler should be disposed immediately after used.




Recommended Storage Condition

- Store at dry, clean, and well-ventilated area at room temperatures between 15° to 25°C (59° - 77°F).
- Do not expose the product to sunlight, UV light, as it will change the properties of the product.
- Keep away from moist and store at dry place with humidity level between 30-60% RH.

Recommended Disposal Method

- If it is regular waste, dispose of the product into a medical pulp macerator, if available. Otherwise, dispose of the product in accordance with the local waste protocol.
- If the product or content is highly contaminated or contagious, treat it as a clinical waste and dispose it in accordance with clinical waste protocol, i.e. incinerate.

Product symbol



Kabomed GmbH
 Freiheitsstraße 43/1
 A-2331 Vösendorf, Austria
office@kabomed.at

UKRP
Curas Ltd.
 Country Gates, Ashton Rd.
 Bristol BS3 2JH, England
info@curas.com



Curas Sdn. Bhd.
 69A-B Raja Uda Business Ctr Rd 4
 12300 Penang, Malaysia
info@curas.com