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Patient Information Leaflet Revamil® Melginate

The Product
Revamil® Melginate –Alginat dressing met honing
5x5 cm / 10x10 cm honey-impregnated alginat sheet for the treatment of wounds

Product description
Revamil® Melginate is a calcium alginat sheet impregnated with pure honey that has been produced under controlled conditions. The carefully controlled production process guarantees that Revamil® Melginate with honey does not contain any traces of pesticides and that the same quality can be reproduced consistently. When Revamil® Melginate is applied to exuding wounds a gel is formed by absorption of wound fluid which protects the wound against maceration. Due to the properties of the alginat applied, i.e. the alginat is fast gelling, non-toxic and has a high wet integrity, the alginat gel covers the wound and maintains a moist wound environment. The alginat sheet impregnated with honey, before and after gelling, protects the wound against bacterial ingress. The gel facilitates easy dressing change because it does not stick to wound and no residues remain in the wound after removal.

Instructions for use
1. When required first clean the wound in the conventional way.
2. Open the outer packaging of the pouch.
3. Wear disposable gloves when handling the alginat sheet.
4. Remove the dressing from the transparent blister.
5. Apply dressing to wound surface.
6. Cover and secure Revamil® Melginate dressing with an appropriate non-occlusive secondary dressing.
7. In case of high exudate level change the dressing every 24 hours.
8. In case of moderate exudate level change the dressing every 3 days or earlier when the dressing has reached its absorbance capacity.
9. If part of the alginat dressing adheres to the wound, moisten the dressing with sterile water or normal saline solution.
10. Gently discard the alginat dressing from the wound and discard.

Indications
Revamil® Melginate impregnated with honey is intended for the management of different types of exuding wounds, including but not limited to, different kinds of ulcers, infected wounds, surgery and radiation induced excoriated wounds, minor (1st and 2nd degree) burns.

Contra indications
• Do not use the product if the patient is known to be allergic or hypersensitive to honey. An allergic reaction or hypersensitivity can cause redness, irritation and possibly minor pain in the skin or wound.
• Do not use on a dry wound. Using Revamil® Melginate on a dry wound could lead to lower performance of the product and slower wound healing.
• In the event of unexpected side effects after use, contact your doctor or the manufacturer. Side effects can cause inconvenience to the patient.

Side effect
During the first minutes after application, the product can cause a light stinging sensation. The stinging sensation is probably caused by the low pH and high conductivity of the honey and has no effect on the performance of the product.

Safety precautions
• Do not use if honey is leaking from the inner packaging. Leakage can lead to contamination of the product with bacteria and to wound infection. Because of the antibacterial protection that honey provides, the chance of wound infection is considered very low.
• Do not use if the packaging (pouch) is not tightly sealed. Packaging that is not tightly sealed can lead to non-sterility of the wound dressing and to wound infection. Because of the high antibacterial protection that honey provides, the chance of wound infection is considered very low.
• Do not apply this product for more than 30 days continuously. Treatment for more than 30 days continuously could lead to inconvenience to the patient, such as irritation and/or minor pain.
• Do not reuse an already used dressing. Reuse of a dressing can lead to loss of integrity of the honey, contamination of the product and/or contamination or infection of the wound.
• Do not reutilize an opened packaging. Reutilization of an opened packaging can lead to loss of function of the honey, contamination of the product and/or contamination or infection of the wound.
• Do not use the product after the expiration date as this can lead to lower performance of the product and to inconvenience for the patient.

Storage conditions
Store between 10 and 30 °C.

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