Odozor Control Dressing / Geschwulstabsorbierendes Wundauflage / Meccazione per il controllo degli odori / Pansement destiné au contrôle des odeurs / Apoio de control do odor / Pense para Control do Odoz / Geurabsorberende Wondauflage / Högtryck Absorberande bandage / Luftabsorbierende bandage / Enfiado Odoroso Odozor / Enfiado Odoroso Odozor

PRODUCT DESCRIPTION
Carborex® Odozor Control Dressing is a sterile fly-free non-adhesive dressing with an absorbent wound contact layer (containing Alginate and Hydrofibre™) an activated charcoal central pad, an absorbent padding layer, and self-flocking outer contained within a peel pouch.

INTENDED USE
Carborex® Odozor Control Dressing is intended for the management of malodorous wounds. It may be used as a primary or secondary dressing as per the indications for use.

INTENDED USERS
Carborex® Odozor Control Dressing is intended to be used by healthcare professionals and by carers and patients under the direction of a Healthcare Professional.

TARGET PATIENT POPULATION
Carborex® Odozor Control Dressing is intended for use in adults and children. It is not recommended for infants less than 12 months of age.

CLINICAL BENEFITS
The clinical benefits of Carborex® Odozor Control Dressing are:
• Management of malodorous wounds
• Absorb and control exudate
• To support maintaining a moist environment for wound healing

INDICATIONS
Carborex® Odozor Control Dressing is indicated for the management of pressure ulcers, venous leg ulcers, arterial leg ulcers, mixed aetiology ulcers and diabetic foot ulcers.

CONTRAINDICATIONS
Carborex® Odozor Control Dressing should not be used on individuals who are sensitive to or who have had an allergic reaction to the dressing or any of its components including any of the components of the adhesive or other seaweed derivatives.

PRECAUTIONS AND OBSERVATIONS
• The use of Carborex® Odozor Control Dressing is not recommended for infants less than 12 months of age due to the presence of manganese within the dressing (120ppm). Rarely, irritation and maceration of the skin and hypopygia around the dressing have been reported.
• Sterility is guaranteed unless the pouch is damaged or opened prior to use. Do not use the device if the packaging is damaged or open prior to use and dispose of the device according to local disposal regulations.
• Do not use after stated expiry date.
• The dressing may be used on infected wounds only under the care of a healthcare professional.
• Carborex® Odozor Control Dressing is not compatible with petroleum-based products.
• Carborex® Odozor Control Dressing is not designed to cut.
• Reseal may lead to increased infection risk. Do not attempt to re-seal.
• Physical properties of the device may change over time for optimal use for intended use.
• Should you observe evidence of maceration, irritation, inflammation (itching of skin, hypergranulation (excess tissue formation) or sensitivity (allergic reaction), contact a healthcare professional.
• After use, this product may be a potential biohazard. Discard as per local clinical practice.
• If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to the competent authority of the Member State in which the use of the device was established.

ASSOCIATED PRODUCTS & THERAPIES
Carborex® Odozor Control Dressing may be used as a primary or secondary dressing if the use of Carborex® Odozor Control Dressing is not recommended for infants less than 12 months of age due to the presence of manganese within the dressing (120ppm). Rarely, irritation and maceration of the skin and hypopygia around the dressing have been reported.

CONTRACEPTION & FERTILITY
Carborex® Odozor Control Dressing can be worn for up to 3 days. Dressings should be changed earlier if clinically indicated. It is safe to continue using Carborex® Odozor Control Dressing for as long as clinical practice allows.

DIRECTIONS FOR USE
Wound site preparation and skin cleansing:
These dressings are sterile and should be used on a clean wound site. If required, the wound site should be debrided and necrotic tissue removed. Once the wound site is clean, well and dry the surrounding skin.

Dressing preparation and application:
Do not cut dressing.
1. Choose a dressing sized so that it is larger than the wound edge leaving a sufficient border for shallow wounds the dressing may be placed directly onto the wound as a primary dressing and for cavity wounds Carborex® can be laid over a wound filler gel as a secondary dressing.

2. Place foam (non-foam) surface over wound or primary dressing.
3. Secure Carborex® in place by pulling the adhesive border over the surrounding skin.
4. The absorbent wound contact layer will take up exudate and form a self-gel.
5. Change the dressing when clinically indicated or if there is exudate visible through the top layer of the dressing or if the dressing is being absorbed. With Carborex® malodorous wounds, Carborex® may be left undisturbed for up to three days.

Removal and disposal:
Carborex® lift away from the wound by grasping one corner of the dressing. Removal from lightly soiled wounds may be assisted by saturating the dressing with sterile normal saline.
2. After use, this product may be a potential biohazard. Discard as per local clinical practice.

STORAGE PRECAUTIONS
Store at room temperature (10-25°C/50-77°F).
Avoid refrigeration and exposure to direct sunlight. Keep Dry.
If further information is required, please contact Convatec Professional Services.
© 2025 Convatec
TM indicates a trademark of the Convatec Group

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