[illegible]

GB	B. Braun Medical Ltd. Brookdale Road Sheffield, S15 2PW United Kingdom	GB	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	GB	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro	GB	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
DE	B. Braun Melsungen AG 34122 Melsungen - Germany	DE	B. Braun Medical A/S 2000 Frederiksberg - Denmark	DE	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	DE	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
ES	B. Braun Medical, S.A.U. Carrereta de Terrassa, 121 08191 Rubí (Barcelona) - España	ES	B. Braun Medical Oy 00350 Helsinki - Finland	ES	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	ES	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
FR	B. BRAUN MEDICAL SAS 92210 Saint-Cloud - France	FR	B. Braun Medical Oy 00350 Helsinki - Finland	FR	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	FR	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
IT	B. Braun Milano SpA. 20161 Milano - Italy	IT	B. Braun Medical Oy 00350 Helsinki - Finland	IT	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	IT	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
NL	B. Braun Medical B.V. Eckerdorp 10, NL-5340 AP Oss. Netherlands	NL	B. Braun Medical Oy 00350 Helsinki - Finland	NL	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	NL	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
NO	B. Braun Australia Pty. Ltd. Level 5, 7-9 Irvine Place Bella Vista NSW 2153 Phone: 1800 251 705	NO	B. Braun Medical Oy 00350 Helsinki - Finland	NO	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	NO	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
PT	Estrela Consiliger Pedrosa 80, Queluz de Baixo, 2730-055 Barcelena Portugal	PT	B. Braun Medical Oy 00350 Helsinki - Finland	PT	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	PT	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
RU	Б.Браун Медикал ФООД г. София 1592 - България Тел. +359 2 807 6740	RU	B. Braun Medical Oy 00350 Helsinki - Finland	RU	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	RU	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
SE	B. Braun Medical AB Box 110 - 182 12 Danderyd - Sweden	SE	B. Braun Medical Oy 00350 Helsinki - Finland	SE	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	SE	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
SI	B. Braun Medical B.V. Eckerdorp 10, NL-5340 AP Oss. Netherlands	SI	B. Braun Medical Oy 00350 Helsinki - Finland	SI	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	SI	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
SK	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	SK	B. Braun Medical Oy 00350 Helsinki - Finland	SK	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	SK	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
SV	B. Braun Medical AB Box 110 - 182 12 Danderyd - Sweden	SV	B. Braun Medical Oy 00350 Helsinki - Finland	SV	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	SV	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
TH	Б.Браун Медикал ФООД г. София 1592 - България Тел. +359 2 807 6740	TH	B. Braun Medical Oy 00350 Helsinki - Finland	TH	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	TH	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
TR	B. Braun Medical B.V. Eckerdorp 10, NL-5340 AP Oss. Netherlands	TR	B. Braun Medical Oy 00350 Helsinki - Finland	TR	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	TR	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro
UK	B. Braun Medical Ltd. Brookdale Road Sheffield, S15 2PW United Kingdom	UK	B. Braun Medical Oy 00350 Helsinki - Finland	UK	B. Braun Medical s.r.o. 148 00 Praha 4 • Česká republika	UK	B. Braun Medical S.R.L. S. Sandares, Str. Berind nr.1 Jed. Timiș. 207375 - România Tel. +40 256 284 905 contact@bbrain.com www.bbrain.ro

5120808 / 2023-10 / N1





B. BRAUN MEDICAL SAS | 26 Rue Armengaud | 92210 Saint-Cloud | FRANCE

Document No.: NA - Version: 2.0 - Document ID: GMS-OMSC-000534 Print Date - Gedruckt am: 2023-12-22 19:15 (CET)


Flexima®/Softima® Active

Effective

Flexima®/Softima® Active

EN Instructions for Use

What Flexima® / Softima® Active is and what it is used for:
Flexima® / Softima® Active is a 1-Piece collection system for enterostomates or urostomates comprising of different types of bags:
■ Drainable for colostomates with solid or liquid stools, ileostomates and jejunostomates with secret and/or High Flow cap (available only for Active O'Connel) allowing connection to a B. Braun distal bag.
■ Drainable for urostomates with a "Jip" tap and ure connector allowing connection to a B. Braun distal bag.
The soft cornex-skin protector Flexima® / Softima® Active O'Connel is specifically designed to ensure a secure seal around stomas and is suitable for patients with shallow to moderately deep skin folds, thin or slightly retracted.
Flexima® / Softima® Active may be used concomitantly with B. Braun stoma care accessories.
Flexima® / Softima® Active device is made of hydrocolloid skin adhesive, polymer/polymer plastics and nonwoven film.
Intended use:
Flexima® / Softima® Active is intended for collection and containment of stools or urine by patients with enterostomy or urostomy.
The product is not intended to be used on injured skin or mucous membrane.
Indication:
The product is indicated for patients with an enterostomy and urostomy.
Flexima® / Softima® Active can only be used in adults. Not for children. Not for infants.
The intended users are healthcare professionals (HCP), caregivers and patients.
The responsible HCP must ensure that the patients and caregivers are duly trained and instructed on the use of ostomy bags, peristomal skin complications and leakage.
Contraindications:
■ Medical contraindications (see product description)
■ Inability to self-stoma care in the absence of an appropriately trained caregiver or attendant for patient with poor manual dexterity and cognitive impairment.
Side effects:
■ Unknown hypersensitivity reaction to any material used leading to anaphylactic or anaphylactoid reactions.

WARNING: If severe hypersensitivity reactions occur, remove the product immediately, discontinue its use and contact a physician.
WARNING: Leakage can promote peristomal skin complications. Therefore use of stoma care accessories to promote skin health can be considered.
1 Preparing the skin
Wash your hands before starting the procedure.
Gently clean the skin around the stoma with lukewarm water. Avoid ether, alcohol or any other topical product which could interfere with the skin protector adhesion. Carefully dab dry the skin with soft and non-shedding cloths.
Choose the appropriate size of the skin protector bag to protect the peristomal skin from exposure to stool and urine. Stoma size may vary from day to day. The skin protector should be adapted accordingly for appropriate fitting of skin protector bag.
2 Preparing the skin protector
If you are using a Cut-to-Fit product, cut the central opening of the hydrocolloid adhesive to the size of your stoma (A) just before the application, possibly with the help of the Stoma Guide. Do not cut beyond the maximum possible diameter and make sure that there is always a 2 to 3 mm margin around your stoma in order not to injure the skin.
Note: Pre-Cut products are not intended to be cut out.
3 Preparing the bag
Separate the bag walls by curling it slightly to introduce air into it. For drainable "Roll Up" bags, close the closure outlet:
■ Drainable "Roll Up" bags: Fold the extremity of the sleeve four times using the guides and fold the tab back (B). To seal the outlet, make sure that the two self-fastening strips are aligned, facing each other.
■ Drainable "High Flow" bags: Once the cap as far as maximum, turning it clockwise. (LX)
4 Applying the skin protector
Peel off the release (B). Position the bottom of the skin protector around the base of the stoma. Press the bottom part of the skin protector smoothly onto the skin, then place and press onto the lower part (C). The skin protector is flexible and adapts perfectly to body contours. (C)
In case of error in the positioning, it is possible to remove the bag

MD	EN	Medical device as per 2017-745/EU	EN	Distributor
DE	DE	Medizinprodukt gemäß 2017-745/EU	DE	Handler
ES	ES	Producto sanitario según el Reglamento (UE) 2017-745	ES	Distribuidor
FR	FR	Dispositif médical selon 2017-745/EU	FR	Distributeur
IT	IT	Dispositivo medico secondo 2017-745/EU	IT	Distributore
NL	NL	Medisch hulpmiddel conform 2017-745/EU	NL	Distributeur
BG	BG	Медицинско изделие съгласно 2017-745/EU	BG	Дистрибутор
CS	CS	Zdravotnický prostředek dle nařízení (EU) 2017/745	CS	Distributor
DA	DA	Medicinsk udstyr i henhold til 2017-745/EU	DA	Distributør
EL	EL	Ιατροτεχνολογικό προϊόν σύμφωνα με τον Κανονισμό 2017/745/ΕΕ	EL	Διανομέας
ET	ET	Meditsiinsesale 2017-745/EU kohast	ET	Levitaja
FI	FI	Lääkinnällinen läite direktiivin 2017-745/EU mukaisesti	FI	Jakelija
HU	HU	Orvosi eszköz, a 2017/745/EU rendeletnek megfelelően	HU	Forgalmazó
JA	JA	2017-745/EUに倣った医療機器	JA	販売代理店
LV	LV	Medicīniskā ierīce atbilstoši ar 2017-745/ES	LV	Implantājs
PL	PL	Wyrob medyczny zgodnie z rozporządzeniem (UE) 2017/745	PL	Dystrybutor
PT	PT	Dispositivo médico de acordo com a diretiva 2017-745/EU	PT	Distribuidor
RO	RO	Dispozitiv medical conform 2017-745/EU	RO	Distributor
SK	SK	Медицинское изделие согласно 2017-745/EU	SK	Дистрибутор
SL	SL	Zdravotniška pomočila podla norme 2017-745/EU	SL	Distributer
SR	SR	Medicinski pripomoček u skladu sa 2017-745/EU	SR	Distributer
SV	SV	Medicinsko sredstvo i skladu sa 2017-745/EU	SV	Distributör

EN Instructions for Use

and replace it properly.

In order to improve the adhesion of Flexima® / Softima® Active O'Convex products, place hands over the skin protector to warm it up.

5 Emptying of drainable bag

Position yourself appropriately.

■ Drainable "Roll'Up" bag:

Release the two self-fastening strips and unroll the sleeve toward the toilet bowl. Be careful, emptying will start right at this point. To fully open the outlet, apply slight pressure with your fingers on guide (E). Once the emptying is complete, expel, as much as possible, any matter still in the sleeve by pressing it lightly from top to bottom. Carefully clean its extremity (use only non-shedding material) (F).

■ Drainable "High Flow" bag:

Remove the cap by pulling it straight and point it towards the toilet bowl. Once the emptying is complete, clean the cap.

■ Drainable "Uro" bag

Point the tap towards the toilet bowl, slowly turn the tap counter-clockwise until the desired flow rate is obtained.

After emptying of drainable bag, close the closure outlet again as indicated above. (J,K)

6 Using B. Braun belt (for Flexima® / Softima® Active O'Convex only)

The soft convex skin protector helps to prevent leakages by tight-fitting the peristomal area. For extra security you can use one out of the three positions of the belt loops of the soft convex skin protector to fix the belt. Choose two loops that are diametrically opposed. (H, I)

7 Using filter protection stickers

Before bathing or showering, stick the filter protection sticker over the gas evacuation openings in order to maintain filter performances (D). For Split cover, this must require lifting the comfort film and to stick the filter protection filter directly on the bag films. After bath or shower, remove the filter protection stickers.

8 Compatibility with urostomy B. Braun distal bags

B. Braun offers a complete solution with urine collection bags for long-time intervals and during the night, such as Uribag® and Urimed® that are compatible to all B. Braun Uro bags.

The connection of the bag to a urostomy B. Braun distal bag can require the use of the supplied uro connector. To clip on the uro connector or the B. Braun Uribag®, press the two wings and insert it into the bag tap. Open the tap turning it half a turn counter-clockwise to allow urine to flow in the B. Braun distal bag. To un-clip, close the tap of the bag by turning it half clockwise, pull out the uro connector or the B. Braun Uribag® by pressing the two wings.

9 Compatibility with Flow Collector®

B. Braun is facilitating the management of high output stoma with B. Braun Flow Collector® which is compatible with all the Flexima® / Softima® range (High Flow products). It will reduce the need to get up at night to empty the system.

10 Removing the bag

Change/empty the bag if volume capacity is reached and/or if you experience odors or ballooning. Diet can affect the stool consistency, ballooning, and odor.

Pull slightly on the top of the skin protector and remove the skin protector from the top to the bottom while maintaining gentle pressure on the skin around the stoma with your other hand. Always remove the skin protector/ bag gently and carefully to reduce peristomal skin complication and additionally to reduce discomfort and pain.

Discard the ostomy bag according to the current local regulations.

Do not flush down the toilet.

Clinical benefits & performances:

Flexima®/ Softima® Active allows the collection and containment of stools or urine.

If you note any change on your stoma and/or peristomal skin, please contact your health care professional.

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