Bard® Flip-Flo™ Catheter Valve Instructions for Use

Description: The Flip-Flo™ Catheter Valve is a device that can be added to a catheter that is placed in the urethra or in suprapubic region. The valve is placed on the end of the catheter and can be manually opened or closed to prevent urine outflow. Sólo para pacientes adultos.

Indication for Use: For use with an indwelling catheter.
Intended Use: The Flip-Flo™ Catheter Valve is intended to facilitate urine drainage and/or collection. For urological use only.
Intended Users: Qualified Healthcare Professionals, Caregiver, Patient/User.

Contraindications:
• Reduced bladder capacity • Cognitive impairment • No bladder sensation
• Inufficient manual dexterity to operate the Flip-Flo™ Catheter Valve

Warnings: This product should not be used without assessment of bladder function or an appropriate medical professional.
This is a single use device. Do not re-sterilize any portion of this device. Reuse and/or repackaging may create risk of patient or user infection. The manufacturer does not warrant the integrity of any material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.

Leave Flip-Flo™ Catheter Valve in open position to allow continuous urine drainage overnight. It is recommended that the Flip-Flo™ Catheter Valve is changed every 5-7 days. Simply disconnect the valve from the catheter and discard.

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

If you have any clinical concerns with the use of this device for your condition, please contact your healthcare professional.

Users and/or patients within the European Union, should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user/patient is established. Users and/or patients of the European Union should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the country in which the user and/or patient is established.

Precautions: Do not use device if package is opened or damaged.
Instructions for Use:
Warning: This product should not be used without assessment of bladder function by an appropriate medical professional.

Instructions:
• Wash hands thoroughly before and after operating Flip-Flo™ valve.
• Do not touch the catheter before or after the valve connector.
• Allow urine to flow freely as advised by your doctor or nurse.
• To open the Flip-Flo™ valve, pull lever down.
• Allow urine to drain into tray or on appropriate receptacle.
• To close the Flip-Flo™ valve, pull lever up.
• A night drainage bag may be attached to the flexible connector to allow continuous urine drainage overnight.

Warning: Leave Flip-Flo™ valve in open position.
It is recommended that the Flip-Flo™ valve is changed every 5-7 days. Simply disconnect the valve from the catheter and discard. Attach as follows following the Use Instructions.

For U.S. patents that may apply, see [BD.com/patents](#).

The BD English symbols glossary is located at [bd.com/en-us/symbols-glossary](#).

Standard Ref. No.
Explanatory Text

Manufacturer 5.1.1 Indicates the medical device manufacturer.

Authorized Representative in the European Community 5.1.2 Indicates the authorized representative in the European Community.

Date of manufacture 5.1.3 Indicates the date when the medical device is fabricated.

Use by 5.1.4 Indicates the date after which the medical device is not to be used.

Lot number 5.1.5 Indicates the manufacturer's batch code so that the batch or lot can be identified.

Catalogue number 5.1.6 Indicates the manufacturer's catalogue number so that the medical device can be identified.

Sterilized using ethylene oxide 5.2.3 Indicates a medical device that has been sterilized using Ethylene Oxide.

Do not sterilize 5.2.6 Indicates a medical device that is not to be sterilized.

Do not use if package is damaged and consult instructions for use 5.2.8 Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional instructions for use.

Single sterile barrier system 5.2.11 Indicates a single sterile barrier system.

Single use 5.4.2 Indicates a medical device that is intended for one use, or for use on single patient during a single procedure.

Consult instructions for use or consult electronic instructions for use 5.4.3 Indicates the need for the user to consult the instructions for use.

Medical device 5.7.7 Indicates the item is a medical device.

Unique device identifier 5.7.10 Indicates a carrier that contains unique device identifier information.

Packaging unit 2794 To indicate the number of pieces in the package.

Open here 3079 To identify the location where the package can be opened and that the method of opening it.

Authorized Representative for Switzerland N/A Indicates the authorized representative for Switzerland.

Authorized Representative for the United Kingdom N/A Indicates the authorized representative for the United Kingdom.

Not Made with Natural Rubber Latex N/A To indicate the device is Not Made with Natural Rubber Latex.

BD, Bard and Flip-Flo are trademarks of Becton, Dickinson and Company or its affiliates. © 2023 BD. All rights reserved.

C. R. Bard, Inc.
8195 Industrial Boulevard
Covington, GA 30014 USA • +1.844.823.5433 [bd.com](#)

Becton Dickinson Ireland Ltd.
Donore Road, Drogheda, Co. Louth, A92 YW26
Ireland

BD Switzerland Srl
Route de Crassier 17, Business Park Terre-Bonne
Bâtiment A4, 1262 Eysins, Suisse

Becton Dickinson UK, LTD
1030 Eskdale Road, Wimshurst Triangle
Wokingham, Berkshire, England, RG41 5TS

Australization Sponsor:
Becton Dickinson Pty Ltd, TIA Bard Australia Pty Ltd
66 Waterloo Road, Macquarie Park, NSW 2113
Australia

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Instructions:
• Wash hands thoroughly before and after handling the Flip-Flo™ valve.
• Do not touch the catheter before or after the valve connector.
• Allow urine to flow freely as advised by your doctor or nurse.
• To open the Flip-Flo™ valve, pull lever down.
• Allow urine to drain into tray or on appropriate receptacle.
• To close the Flip-Flo™ valve, pull lever up.
• A night drainage bag may be attached to the flexible connector to allow continuous urine drainage overnight.

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