DAR™

Breathing System

en	Instructions For Use
fr	Mode d'emploi
it	Istruzioni per l'uso
de	Gebrauchsanweisung
es	Instrucciones de uso
sv	Bruksanvisning
nl	Gebruiksaanwijzing
pt	Instruções de uso
da	Brugsanvisning
fi	Käyttöohjeet
el	Οδηγίες χρήσης
cs	Návod k použití
hu	Használati útmutató
ru	Инструкция по применению
pl	Instrukcja użytkowania
tr	Kullanma Talimatı
no	Bruksanvisning

		Rated flow	Resistance to flow
ADT	22mm	60 l•min ⁻¹	≤0,2 kPa (2 cm H ₂ O)
PED	15mm	30 l•min ⁻¹	≤0,2 kPa (2 cm H ₂ O)
INF	10mm	15 l•min ⁻¹	≤0,2 kPa (2 cm H ₂ O)


Fig. 1


Fig. 2

Fig. 3


Fig. 4

CLEAN/STERILE (CHECK LABEL) AND READY FOR USE: unless package has been opened or damaged. Single use device to be used on a single patient. Discard after use.

Identification of a substance that is contained or present within the product or packaging.

Identification of a substance that is not contained or present within the product or packaging.

This product contains DEHP. When used as indicated, very limited exposure to trace amounts of DEHP may occur. There is no clear clinical evidence that this degree of exposure increases clinical risk. However, in order to minimize risk of DEHP exposure in children and nursing or pregnant women, this product should only be used as directed. Refer to primary labeling.

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

DESCRIPTION:

Breathing circuit for anaesthesia (Fig. 1), intensive care (Fig. 2) and/or CPAP support/breathing therapy (Fig. 3) for adult, paediatric and neonatal patients (indication on label). The circuit can be provided with additional ports for the monitoring of temperature and pressure (points A and B Fig. 2) and one-way valves (i.e. PEEP, Fig. 4). An anaesthesia mask with inflatable cushion can be also connected to the breathing circuit (point A, Fig. 1).

INDICATIONS:

The circuits are intended to connect the patient to a breathing gas source during anaesthesia or intensive care.

The devices are intended for single use in adult, paediatric and neonatal patients undergoing mechanical ventilation in anaesthesia and intensive care.

The circuits are in compliance with ISO 5367 and EN 12342 standards requirements.

The conformity evidence of the circuits to these standards, concerning resistance to flow values, is in the above chart.

CONTRAINDICATIONS:

Do not use the devices in positions different from those described in the picture.

DIRECTIONS FOR USE:

1. Remove the circuit from its protective package.
2. Each circuit component must be visually inspected and functionality and patency must be carefully checked immediately before its use on patient.
3. Check circuit integrity when under pressure.
4. Connect the circuit to the ventilator and to the patient.

Pay particular attention to the following indications:

- a. in presence of a pressure monitoring line, this must be connected to the pressure port on the patient Y connector and to the pressure port on the ventilator. If the line is not required, use the appropriate plug to seal the pressure port on the patient Y connector.
- b. If a temperature monitoring is required, remove the plug from the port and install the temperature control probe. Seal the port with its plug when the temperature probe is not used.
- c. CPAP circuits (at continuous flow): in the patient three-ways adapter (T or Y scheme), the expiratory port, opposite to the feeding flow, must always be kept free; its occlusion or partial reduction can induce serious damages to the patient up to causing his/her death.
- d. If the circuit is provided with a water-trap, position it at the lowest point of the circuit limb (inspiratory or expiratory). Monitor the fluid level in the water-trap frequently and empty the device as necessary; make sure the water-trap is properly resealed. Dispose the fluids according to hospital policy.
- e. The circuit can contain one-way valves; verify the proper functionality and flows direction before its use.

- f. If the circuit contains a PEEP valve on the expiratory line, verify its functioning by rotating the regulation ring. The pressure limiting is obtained by rotating anticlockwise the regulation ring. Verify the selected pressure value (from 5 to 20 cm H₂O) using a gauge.
- g. The inflatable mask (point A Fig. 1) is provided with a cushion regulating pressure valve. If necessary, connect the syringe to the valve luer adapter and regulate the cushion inflation up to the desired inflating level. Disconnect the syringe and check the system seal.

Max circuit use: 72 hours. If the breathing circuit contains a filter, filter/HME or HME, replace it every 24 hours. Do not reuse. Discard after use.

Precautions should be taken when discarding the device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

WARNINGS:

- **Prescription device is to be used exclusively under medical supervision.**
- **The product is guaranteed only if the package has not been damaged. Check pack integrity: any damage or opening in the package may compromise sterility and/or performance of the device. In this case, do not use the circuit concerned.**
- **Check the circuit periodically in order to monitor a possible increase of resistance due to condensation into the device.**
- **Performance and safety of the circuit are not guaranteed if used with other accessories.**
- **Ensure that all the connections are correct and firmly secured.**
- **Do not subject the circuit to yield, crushing, or other mechanical stress.**
- **Replace the device whenever deteriorations or anomalies appear (i.e. colour changes, deformations, components detachment, etc..)**
- **In case secretions are detected inside the device, replace it immediately.**
- **Ensure that the circuit is positioned in such a way that it does not form more than one loop or reduce the tube patency.**
- **For any handling of contaminated single use circuits suitable means of protection must be used.**
- **In case a water-trap is used, the overflow of the container could cause a liquid outpouring within the circuit, causing an increase in the resistance to flow potentially dangerous for the patient.**
- **It is not advisable to use a breathing circuit not provided with heated wire, and without a water-trap, along with active humidifiers.**
- **In case an HME, filter/HME or filter is placed on the patient's connector, avoid the concomitant use of an active humidifier and/or a nebulizer positioning upstream the HME, filter/HME or filter since dangerous raisings of the resistance to flow may occur.**
- **To reduce the risk of damage and obstruction by small objects, keep wrapped until required for actual use.**
- **Avoid exposure to elevated temperatures and ultraviolet light during storage of breathing bags both stand alone and with circuits configuration.**
- **Do not use with magnetic resonance imaging (MRI) devices and in any magnetic resonance (MR) environment. Please refer to primary label.**

PRECAUTIONS:

Do not try to refurbish. Do not reuse. Do not soak, rinse, wash, sterilize or treat with any disinfectants (in particular phenolic and alcoholic-based solutions are to be avoided).

POTENTIAL ADVERSE REACTIONS:

Ventilated patients are particularly susceptible to airways infections. Single use devices if not replaced so frequently as indicated may become a potential infections source. Wash accurately the hands before and after the device replacement.

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STERILE EO

Caution: This product contains natural rubber latex which may cause allergic reactions.

Single use

Rx
ONLY

Do not use if package is opened or damaged

Caution, consult accompanying documents

MR unsafe

0123

Part No. 10051316 Rev E 07/2013

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.

© 2011 Covidien.

 Covidien LLC, 15 Hampshire Street, Mansfield, MA 02048 USA.

 Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.

www.covidien.com