Shiley™

Accessories

REF 4SIC	EN US	Instructions For Use
REF 6SIC	FR	Mode d'emploi
REF 8SIC	DE	Gebrauchsanweisung
REF 10SIC	IT	Istruzioni per l'uso
REF 4DCP	ES	Instrucciones de uso
REF 6DCP	PT	Instruções de uso
REF 8DCP	NL	Gebruiksaanwijzing
REF 10DCP	SV	Bruksanvisning
REF CAP	DA	Brugsanvisning
REF DDCP	NO	Bruksanvisning
	FI	Käyttöohjeet
	EL	Οδηγίες χρήσης

This instruction insert applies to the following Shiley Tracheostomy products:

SIC	Spare Inner Cannula
CAP	Cap 15MM
DCP	Decannulation Plug
DDCP	Decannulation Plug Disposable

Identification of a substance that is not contained or present within the product or packaging.

Spare Inner
Cannula

15mm Cap

Decannulation
Plug

Disposable
Decannulation Plug

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

Instructions for Use

POST FOR USE BY ALL TRACHEOSTOMY CARE PERSONNEL

WARNINGS:

CAREFULLY READ THESE INSTRUCTIONS PRIOR TO USE. THESE INSTRUCTIONS ARE APPLICABLE ONLY TO THE SHILEY ACCESSORY PRODUCTS FOR USE WITH THE SHILEY TRACHEOSTOMY TUBES.

STERILE ONLY IF THE PROTECTIVE BAG IS NOT OPEN, DAMAGED OR BROKEN. DO NOT RESTERILIZE. DO NOT EXPOSE TO TEMPERATURES ABOVE 120°F (49°C).

WHEN USING A LASER IN CONJUNCTION WITH THESE TRACHEOSTOMY TUBE ACCESSORIES, CARE MUST BE TAKEN TO AVOID CONTACT OF THE LASER BEAM. SUCH CONTACT, ESPECIALLY IN THE PRESENCE OF OXYGEN-ENRICHED MIXTURES, COULD RESULT IN RAPID COMBUSTION OF THE TUBE WITH HARMFUL THERMAL EFFECTS AND WITH EMISSION OF CORROSIVE AND TOXIC COMBUSTION PRODUCTS, INCLUDING HYDROCHLORIC ACID (HCl).

STERILE. Sterilized with ethylene oxide. For single patient use only.

CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

CAUTIONS:

The Shiley Tracheostomy Tube and Accessory Products are composed of medical grade plastics which do not contain latex material.

Precautions should be taken when discarding any of the tracheostomy tube devices and accessories. Disposal of medical devices should be made in accordance with applicable national regulations for biologically hazardous waste.

The Shiley Tracheostomy Tube is classified as a disposable medical device. Frequent and routine changes of the Tracheostomy Tube are recommended. The manufacturer recommends that the Tracheostomy Tube usage not exceed twenty-nine (29) days. Frequent and routine changes of the Tracheostomy Tube and accessories are recommended and should be evaluated by the attending physician.

Patients in the home care environment should be carefully instructed by a home health care provider in the proper use and handling of the tracheostomy tube and accessory products.

Contents, Descriptions and Intended Use of Shiley Accessory Packages:

SIC - Spare Inner Cannula

Package Contents - ten individually packages units. The SIC is available in four sizes, 4, 6, 8 and 10. The SIC features an integral red 15mm twist lock connector and is designed for temporary use during cleaning or removal of the inner cannula for Shiley Tracheostomy Tubes with twist-lock connectors. Refer to Panel 1 for model types. The Shiley Tracheostomy Tube and the accessory products are intended for use in providing tracheal access for airway management.

SIC SIZE	 (SIC)	 (Outer Cannula)
4	5.0 mm	9.4 mm
6	6.4 mm	10.8 mm
8	7.6 mm	12.2 mm
10	8.9 mm	13.8 mm

Note: The sizing dimensions refer to the I.D. as the inside diameter of the inner cannula at the smallest location. O.D. is the outside diameter of the outer cannula (not including cuff features). Length is the distance from the neck plate to the distal tip on the tube's centerline.

WARNING:
 THE SPARE INNER CANNULA (SIC) WITH RED TWIST-LOCK 15mm CONNECTOR PROVIDED IN THIS PACKAGE CAN ONLY BE USED WITH A SHILEY TRACHEOSTOMY TUBE OF THE SAME SIZE NUMBER.
 DURING AND AFTER ATTACHMENT OF RESPIRATORY OR ANESTHESIA TUBING AND/OR CONNECTORS TO THE SPARE INNER CANNULA, AVOID APPLICATION OF EXCESSIVE ROTATIONAL, LINEAR, OR ROCKING FORCES ON THE TUBING AND/OR CONNECTORS TO PREVENT ACCIDENTAL DISCONNECTION OF THE SPARE INNER CANNULA, OR DAMAGE TO THE TRACHEOSTOMY TUBE.

Instructions

1. Remove the reusable inner cannula and insert the Spare Inner Cannula (SIC) with red 15mm connector.
2. Lock the SIC in place by holding the swivel neck plate firmly with your fingertips and twist the connector clockwise, one-quarter turn past the locking detent.

CAUTION:
To avoid applying pressure against the patient, the neck plate may be stabilized with your free hand during the locking procedure.

Verify that the red twist-lock connector engages securely after each use. If parts become worn or loose, immediately report this to your physician for prompt replacement of the tracheostomy tube.

The SIC is shorter than the original fitted inner cannula and was designed for temporary use. Ten (10) minutes is suggested as the time limit for continual usage. Secretions may build up on the inside of the outer cannula if used longer than the recommended time.

3. The patient's respirations and vital signs should be routinely evaluated after locking the spare inner cannula in the tube.

Cleaning

WARNINGS

- DO NOT USE SOLUTIONS OR CHEMICAL AGENTS OTHER THAN THOSE RECOMMENDED IN THE TABLE BELOW TO CLEAN THE SIC, AS THIS MAY RESULT IN TUBE DAMAGE.
- DO NOT SOAK ANY PART OF THE TUBE IN HYDROGEN PEROXIDE OR ANY OTHER SOLUTION.

• SPARE INNER CANNULA (SIC)	1. CLEAN WITH HYDROGEN PEROXIDE (HALF-STRENGTH), STERILE NORMAL SALINE, OR WATER AND MILD DETERGENT. 2. AFTER CLEANING THOROUGHLY RINSE WITH STERILE SALINE TO REMOVE ANY CLEANING SOLUTION RESIDUES.
-----------------------------	--

Contents, Descriptions and Intended Use of Shiley Accessory Packages:

CAP - 15MM cap

Package Contents - ten individually packaged units. The white 15mm Cap (CAP) is a universal size and is intended to occlude the standard male 15mm connector on all sizes of the Shiley Tracheostomy Tubes. Refer to Panel 1 for model types. The CAP will block all expired and inspired air through the tube, forcing breathing through the mouth and nose, allowing speech. The Shiley Tracheostomy Tube and the accessory products are intended for use in providing tracheal access for airway management.

WARNINGS:

USE THE 15mm CAP ONLY WITH THE SHILEY CUFFLESS FENESTRATED, CUFFLESS, FENESTRATED LOW PRESSURE CUFFED, DISPOSABLE CANNULA CUFFLESS, DISPOSABLE CANNULA CUFFLESS FENESTRATED, AS WELL AS THE SHILEY PEDIATRIC AND NEONATAL TRACHEOSTOMY TUBES.

PRIOR TO ATTACHING THE 15mm CAP MAKE SURE THAT THE FENESTRATION(S) ARE NOT OCCLUDED, THE CUFF IS COMPLETELY DEFLATED AND THAT THERE IS SUFFICIENT AIRWAY FOR THE PATIENT.

1. Establish the patency of the patient's upper airway tract. The patient's airway should be cleared by coughing and/or suctioning before capping the tracheostomy tube with the 15mm Cap.
2. Completely deflate the tracheostomy tube cuff (FEN only).
3. Attach the 15mm Cap on the 15mm tracheostomy tube connector

CAUTION:

To avoid applying pressure against the patient, the neck plate may be stabilized with your free hand during the capping procedure.

Verify that the 15mm Cap is attached securely after each use. If parts become worn or loose, immediately report this to your physician for prompt replacement of the tracheostomy tube.

4. The patient's respirations and vital signs should be evaluated after occluding the tube with the 15mm Cap.

WARNING:

REMOVE THE 15mm CAP IMMEDIATELY IF THE PATIENT HAS DIFFICULTY BREATHING.

5. If the patient needs to be mechanically ventilated while using the 15mm Cap, remove the 15mm Cap, and connected the tracheostomy tube to the ventilator.

or

6. If using the 15mm Cap with the fenestrated inner cannula with green 15mm twist-lock connector, remove the 15mm Cap and fenestrated inner cannula, replace with the non-fenestrated inner cannula with white 15mm connector, inflate the cuff (FEN) and connect the tracheostomy tube to the ventilator.

Cleaning

WARNINGS:

- **DO NOT USE SOLUTIONS OR CHEMICAL AGENTS OTHER THAN THOSE RECOMMENDED IN THE TABLE BELOW TO CLEAN THE 15mm CAP, AS THIS MAY RESULT IN DAMAGE.**
- **DO NOT SOAK ANY PART OF THE TUBE IN HYDROGEN PEROXIDE OR ANY OTHER SOLUTION.**

• 15mm CAP	1. CLEAN WITH HYDROGEN PEROXIDE (HALF-STRENGTH), STERILE NORMAL SALINE, OR WATER AND MILD DETERGENT. 2. AFTER CLEANING THOROUGHLY RINSE WITH STERILE SALINE TO REMOVE ANY CLEANING SOLUTION RESIDUES.
--	---

Contents, Descriptions and Intended Use of Shiley Accessory Packages:

Decannulation Plug (DCP)

Package Contents - ten individually packaged units. The DCP is red for easy identification and is available in four sizes, 4, 6, 8 and 10. The DCP, with twist-lock connector is used to occlude the proximal end of the outer cannula on the Shiley Tracheostomy Tubes. Refer to Panel 1 for model types. The DCP blocks airflow through the tube and directs breathing through the mouth and nose. The Shiley Tracheostomy Tube and the accessory products are intended for use in providing tracheal access for airway management.

WARNINGS:

USE THE DECANNULATION PLUG (DCP) ONLY WITH THE SAME SIZE SHILEY FENESTRATED LOW PRESSURE CUFFED (FEN), AND CUFFLESS FENESTRATED (CFN) TRACHEOSTOMY TUBES.

PRIOR TO INSERTING THE DECANNULATION PLUG (DCP) MAKE SURE THAT THE FENESTRATION IS NOT OCCLUDED, THE CUFF (FEN) IS COMPLETELY DEFLATED AND THAT THERE IS SUFFICIENT AIRWAY FOR THE PATIENT.

1. Establish the patency of the patient's upper airway tract. The patient's airway should be cleared by coughing and/or suctioning before plugging the tracheostomy tube with the Decannulation Plug.
2. Completely deflate the tracheostomy tube (FEN only).
3. Remove the inner cannula with 15mm twist-lock connector and insert the Decannulation Plug with twist-lock connector to the proximal end of the outer cannula.
4. Lock the Decannulation Plug in place by holding the swivel neck plate firmly with your fingertips and twisting the Decannulation Plug clockwise past the locking detent, one quarter turn.

CAUTIONS:

Locking the Decannulation Plug on an angle could strip the locking mechanism.

To avoid applying pressure against the patient, the neck plate may be stabilized with your free hand during the locking procedure.

Verify that the twist-lock connector engages securely after each use. If parts become worn or loose, immediately report this to your physician for prompt replacement of the tracheostomy tube.

5. The patient's respirations and vital signs should be evaluated after plugging the tube with the Decannulation Plug.

WARNING:

IF PATIENT HAS DIFFICULTY BREATHING, IMMEDIATELY REMOVE THE DECANNULATION PLUG, INSERT THE INNER CANNULA AND VERIFY AIRWAY.

6. If the patient needs to be mechanically ventilated while using the Decannulation Plug, remove the Decannulation Plug, insert the non-fenestrated inner cannula with white 15mm twist-lock connector, inflate the cuff (FEN), and connect the tracheostomy tube to the ventilator

Cleaning**WARNINGS:**

- **DO NOT USE SOLUTIONS OR CHEMICAL AGENTS OTHER THAN THOSE RECOMMENDED IN THE TABLE BELOW TO CLEAN THE DCP, AS THIS MAY RESULT IN DAMAGE.**
- **DO NOT SOAK ANY PART OF THE DCP IN HYDROGEN PEROXIDE OR ANY OTHER SOLUTION.**

• DECANNULATION PLUG (DCP)	1. CLEAN WITH HYDROGEN PEROXIDE (HALF-STRENGTH), STERILE NORMAL SALINE, OR WATER AND MILD DETERGENT. 2. AFTER CLEANING THOROUGHLY RINSE WITH STERILE SALINE TO REMOVE ANY CLEANING SOLUTION RESIDUES.
--	--

Contents, Descriptions and intended Use of Shiley Accessory Packages:**Disposable Decannulation Plug (DDCP)**

Package Contents - ten individually packaged units. The DDCP is a universal size and fits any size DCFS, DCFN, DFEN tracheostomy tube. The red DDCP snap-lock connector is used to occlude the proximal end of the outer cannula on the Shiley Tracheostomy Tubes. Refer to Table 1 for model types. The DDCP blocks airflow through the tube and directs breathing through the mouth and nose. The Shiley Tracheostomy Tube and the accessory products are intended for use in providing tracheal access for airway management.

WARNINGS:

USE THE DISPOSABLE DECANNULATION PLUG ONLY WITH THE SHILEY DISPOSABLE FENESTRATED LOW PRESSURE CUFFED, CUFFLESS, AND CUFFLESS FENESTRATED TRACHEOSTOMY TUBES.

PRIOR TO INSERTING THE DISPOSABLE DECANNULATION PLUG (DDCP) MAKE SURE THAT THE FENESTRATIONS ARE NOT OCCLUDED, THE CUFF IS COMPLETELY DEFLATED AND THAT THERE IS SUFFICIENT AIRWAY FOR THE PATIENT.

1. Establish the patency of the patient's upper airway tract. The patient's airway should be cleared by coughing and/or suctioning before plugging the tracheostomy tube with the Disposable Decannulation Plug (DDCP).
2. Completely deflate the tracheostomy tube cuff (DFEN only).
3. Remove the Disposable Inner Cannula (DIC) with 15mm snap-lock connector from the tracheostomy tube.
4. To insert the DDCP, grasp the plug by the pinch clamps of the snap-lock connector and pinch several times to loosen.
5. Holding the pinch clamps, secure one side of the DDCP to the proximal end of the outer cannula.

CAUTION:

Do not attempt to secure both clamps simultaneously. To avoid applying pressure against the patient, the neck plate may be stabilized with your free hand during the locking procedure.

Verify that the Disposable Decannulation Plug (DDCP) is attached securely after each use. If parts become worn or loose, immediately report this to your physician for prompt replacement of the tracheostomy tube.

6. Holding the first side securely, roll the disposable decannulation plug across the connector until the other clamp is in position.
7. Gently but firmly, press the second clamp and let it snap-lock into place. Check for a secure connection.

WARNING:

IF PATIENT HAS DIFFICULTY BREATHING, IMMEDIATELY REMOVE THE DISPOSABLE DECANNULATION PLUG, INSERT THE DISPOSABLE INNER CANNULA AND VERIFY AIRWAY.

8. If the patient needs to be mechanically ventilated while the disposable decannulation plug is in use, remove the disposable decannulation plug, insert the Disposable Inner Cannula (DIC) with 15mm snap-lock connector and then connect it to the ventilator.

CAUTION

The Disposable Decannulation Plug (DDCP) is designed for single use and should not be cleaned or reused.

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© 2011 Covidien.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and/or internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
Covidien Inc, 15 Hampshire Street,
EC REP Covidien Ireland Limited,