

Cystodrain integral Set

Suprapubic drainage sets with Silicone Integrated balloon catheter.

Punction set	Light punction set	Dilation set	Replacement set
Silicone Integral balloon catheter	Silicone Integral balloon catheter	Silicone Integral balloon catheter	Silicone Integral balloon catheter
Splittable trocar (12 cm)	Splittable trocar (12cm)	Dilator	Adhesive tape
Locking sleeve		PTFE Guidewire (90 cm)	Spigot
Adhesive tape		Adhesive tape	Optional PTFE guidewire (90 cm)
Spigot		Spigot	
Scalpel			
Box of 5 - Sterile	Box of 5 - Sterile	Box of 5 - Sterile	Box of 1 or 5 - Sterile

Silicone Integral balloon catheter: Channels: 2 – Eyes: 2 – Balloon size: 3 (Ch.10) / 5 ml – PTFE Guidewire: 0.038" (0.97 mm)



Trocar



Scalpel




Sleeve



Spigot



Dilator



Adhesive tape

Punction sets references

Trocar Ø (Ch/Fr)	14	15	17
Catheter Ø Ch/Fr	10	12	14
Punction set REF.	CD1010	CD1012	CD1014
Light punction set REF.	CD1210	CD1212	CD1214

Dilatation set references

Catheter Ø Ch/Fr	12	14
Set REF.	CD1112	CD1114

Simple replacement sets

Catheter Ø Ch/Fr	12	14	16	18	20
Replacement set REF. (with guidewire)	CD1A12	CD1A14	CD1A16		
Replacement set REF. (without guidewire)	CD1B12	CD1B14	CD1B16	CD1B18	CD1B20
	CD1C12	CD1C14	CD1C16	CD1C18	CD1C20

Box of 1
Box of 5