

PRODUCT INFORMATION* SC® Value PF

Doc 7 Revision:04 Effective: 01/08/2020

Page 1/1

EUROPE						
Product Information:						
Trade Name	SC [®] Value PF					
Туре	Medical Examination Glo	Medical Examination Glove, Polyvinyl Chloride (PVC), Non-sterile, Single Glove, Powder free, Ambidextrous				
Intended Use	To conduct medical examination, diagnostic and therapeutic procedures to protect patient and user from					
	cross contamination or infection.					
Product Conformance	In compliance with European Medical Device Regulation 2017/745. In accordance to: EN 420, EN 455					
Material	Synthetic Poly Vinyl Chlo	Synthetic Poly Vinyl Chloride Paste				
Accelerators	None associated with NRL					
Donning Lubricant	Polyurethane Coating					
Cuff Finishing	Beaded					
Color	Beige					
Design	Ambidextrous					
Storing conditions	Keep dry, Keep away from sunlight					
Temperature limitation	Between 5°C and 25°C					
Residual Powder	No powder lubricant added					
Dimensions of the innerboxes	220 x 115 x 60 mm					
Dimensions of the outer carton	317 x 239 x 233 mm					
Product Specification:	•					
Conform EN 420, EN455						
Reference codes	Size	Order code	Packing unit	Innerboxes		
	Extra Small	02017	1000	10 x 100		
	Small	02012	1000	10 x 100		
	Medium	02013	1000	10 x 100		
	Large	02014	1000	10 x 100		
	Extra Large	02015	1000	10 x 100		
Dimensions	Size	Palm Width (mm)	Length (mm)			
	Extra Small	≤ 80	min. 240			
	Small	80 ± 10	min. 240			
	Medium	95 ± 10	min. 240			
	Large	110 ± 10	min. 240			
	Extra Large	≥ 110	min. 240			
Thickness	Single wall thickness (mm) / Location					
	Palm		0.05			
	Finger		0.05			
Physical Properties	Before Aging	Properties				
	Force at Break (N)	Minimum 3,6 Newton				
	Elongation (%)	Minimum 350%				
Physical Properties	After Aging	Properties				
	Force at Break (N)	Minimum 3,6 Newton				
	Elongation (%)	Minimum 350%				
Shelf life	5 Years upon manufactu	uring date		•		
Performance requirements	Criteria	Inspection Level	AQL			
	1000 ml Water Leak	G-1	1.5			
	Physical properties	S-2	4.0			
	Dimension	S-2	4.0			
	Visual Inspection	S-2	4.0			

Manufactured under ISO9001 Quality Systems

Medica Europe B	V
-----------------	---

Quality assurance Medica Europe operates with a quality management system which complies with the requirements of ISO13485:2016 and the environmental management system ISO 14001:2015

SC® Value and Medica Europe are registered trademarks of Medica Europe B.V.

© Medica Europe. Not to be reproduced or distributed without prior written permission of Medica Europe B.V.

* The product information provided is a guideline of typical performance characteristics of the product and is not to be used as actual product specification.