

Product Data Sheet

Protective clothing



Page 1 (2)

evercare® Essentia isolation gowns REF 2727-01 and 2729-01

Version: 04.07.2016, rev.1

Printed 4.7.2016

| REF | Size | | |
|---|------------------|--|--|
| 2727-01 2729-01 | 120 cm 120 cm | | |
| Colour 2727-01 white, with elasticated cuff 2729-01 yellow, with elasticated cuff [NATURAL LATEX FREE] | | | |
| | | 2727-01 | 2729-01 |
| Raw materials: Body Elastic band at wrists Neck ribbon & belt Sewing thread | | Repellent multilayer SMS PP nonwoven, white or yellow Latex-free PP nonwoven, white or yellow Terylene, white | |
| REACH regulations: | | The products covered by this data sheet specification do not contain carcinogenic, mutagenic or toxic substances, including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA). | |
| Features: | | Soft disposable nonwoven gown with neck ties, waist band and band at wrists. Generous design. | |
| Pack count: | | 2727-01 | 10 pcs/bag, 10 bags (100 pcs)/transport carton |
| | | 2729-01 | 10 pcs/bag, 10 bags (100 pcs)/transport carton |
| Sales unit: | | Transport carton (10 bags of 10 pcs) | |
| Dispenser pouch | | 2727-01 | Plastic bag 330 x 440 mm |
| | | 2729-01 | Plastic bag 330 x 440 mm |
| Transport carton: | | 2727-01 | Double wall corrugated fiberboard 390x295x320 mm |
| | | 2729-01 | Double wall corrugated fiberboard 390x295x320 mm |
| Compliance with legal requirements: | | Medical Devices Directive 93/42/EEC, Class I | |
| Compliance with product standards: | | N/A | |
| Quality standards: | | The product is manufactured in compliance with ISO 13485 | |
| Label information: | | in accordance with EN 15223-1 and EN 1041 | |

Product Data Sheet

Protective clothing







Page 2 (2)

evercare® Essentia isolation gowns REF 2727-01 and 2729-01

Version: 04.07.2016, rev.1

Printed 4.7.2016

| | | |
|--------------------------------|--|---|
| Storage recommendation: | Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight and other intensive light sources. | |
| Disposal recommendation | Follow local instructions | |
| Risk if reused: | The product is single use. If reused performance of product may deteriorate, cross contamination may occur | |
| Barcode / Unit | Dispenser | Transport carton |
| N/A | 6438129321708 2727-01 dispenser  6 438129 321708 | 6438129521702 2727-01 carton  6 438129 521702 |
| N/A | 6438129321715 2729-01 dispenser  6 438129 321715 | 6438129521719 2729-01 carton  6 438129 521719 |