

UROMED »SUPRASOFT« DOPPEL-BALLON-KATHETER

REF 2960

Ohne Spitze, zentral offen
1 Auge, Länge 36 cm

Doppel-Ballon-Katheter zur suprapubischen Blasendrainage

Der vordere Ballon des von Prof. Walz entwickelten Katheters verhindert durch die große Auflagefläche Irritationen an der Blasenwand und bietet sich damit für den Einsatz bei sensiblen Patienten an. Die Drainage erfolgt durch das Auge zwischen den Ballons.

- besonders geeignet bei Tenesmen
- Blockung beider Ballons über einen Kanal
- großlumig
- Röntgenkontraststreifen
- Soft-Ballons
- reines Silikon
- halbfest

Double Balloon Catheter for suprapubic drainage of the urinary bladder

The catheter, developed by Prof. Walz, with a large contact surface of the distal balloon prevents irritations of the bladder wall and is destined for sensible patients. The drainage will be effected through the eye between the balloons.

- without tip, central open, 1 eye, length 36 cm
- especially for tenesmen
- both balloons to be blocked by using one channel
- wide lumen
- x-ray contrast line
- soft balloons
- pure silicone
- semi-solid



REF REF	Größen Sizes	Ballon Balloon	PZN PZN
2960-14	Ch. 14	15-20 ml	04741880
2960-18	Ch. 18	15-20 ml	04741897

Verpackung: doppelt, steril
Original-Karton: 10 Stück

Packing: double, sterile
Original carton: 10 pieces

UROMED
PRODUKTE FÜR DIE UROLOGIE