

PRODUCT

LEG URINE BAG WITH NON-RETURN VALVE AND PUSH PULL OUTLET

Bag's body:	made of a-toxic, soft and strong PVC DEHP FREE
Bag's color:	transparent
Bag's Tubing:	made of a-toxic, transparent PVC DEHP FREE PVC granule's high quality gives a high «memory» effect to the tubing making it very folding enduring. The tubing length could be adapted according to the customers need (between 10 cm and 50 cm).
Connector:	universal, white, made of a-toxic material (ABS);
Non Return Valve:	close, to prevent urine and gas back-flow
Tap and emptying:	PUSH PULL outlet tap; it allows the bags to be easily and quickly emptied with the guarantee of preventing obstructions during the outflow
Welding:	double, excellent holding.
Attachment point:	2, strengthened to guarantee a solid and stable joining, by the fastening system, to the bed and to assure the maximum security in the usage and movements of the bag.
Graduation:	asymmetrical to allow a correct reading of the liquid in
Bag's Capacity:	750 ml or 500 ml
Trademark:	Securdrain Vià!
Sterilization:	by ETO+CO2 blend
Custom Tariff:	39232910

Codes and Package

CODE	DESCRIPTION	PACKAGE
001019133350	Leg Bag, 750 ml, tubing 35 cm, with non return valve and Push Pull outlet. BULK, not sterile	400 pcs / carton Carton Dimensions: 0.39*0.29*0.39 cm
001013133350	Leg Bag, 750 ml, tubing 35 cm, with non return valve and Push Pull outlet. Package of 10 pcs – equipped with 2 straps each polybag	40 polybag of 10 pcs per carton. 400 pcs / carton Carton Dimensions: 0.39*0.29*0.39 cm
001013133100	Leg Bag, 750 ml, tubing 10 cm, with non return valve and Push Pull outlet. Package of 10 pcs – equipped with 2 straps each polybag	40 polybag of 10 pcs per carton. 400 pcs / carton Carton Dimensions: 0.39*0.29*0.39 cm
001014133350	Leg Bag, 750 ml, tubing 35 cm, with non return valve and Push Pull outlet. Package of 30 pcs – equipped with 2 straps each polybag	12 polybag of 30 pcs per carton. 360 pcs / carton Carton Dimensions: 0.39*0.29*0.39 cm
001014133100	Leg Bag, 750 ml, tubing 10 cm, with non return valve and Push Pull outlet. Package of 30 pcs – equipped with 2 straps each polybag	12 polybag of 30 pcs per carton. 360 pcs / carton Carton Dimensions: 0.39*0.29*0.39 cm
001013134351	Leg Bag, 500 ml, tubing 35 cm, with non return valve and Push Pull outlet. Package of 10 pcs – equipped with 2 straps each polybag	40 polybag of 10 pcs per carton. 400 pcs / carton Carton Dimensions: 0.39*0.29*0.39 cm

Technical Drawing

