



Information

Silicone backed comfort fabric straps can also be ordered separately (see Accessories for catheterisation*).



Three-chamber leg urine collection bags

Discretion and comfort

The Three-Chamber Leg Urine Collection Bag compared to the one-chamber leg bag is designed to give the mobile user even more discretion.

Inconspicuous

Due to the three interconnected chambers, the collector:

- Adapts to the leg
- Can fill up evenly
- Is worn discreetly
- Is hardly visible underneath clothing

Triple convenience

- The ridged universal connector with protective cap helps to enable a secure connection with the catheter funnel.
- The needle-free urine sampling port* designed to allow urine withdrawal with a Luer-slip syringe without the risk of needle-stick injury
- The non-reflux valve is designed to prevent the urine from flowing back.

Triple comfort

- Fleece material on the bag's reverse side
- The non-drip lever tap outlet valve can be operated with only one hand, even with limited manual dexterity – it includes an extension tube for further urine collection bags
- Two silicone backed comfort fabric straps** are included in each packaging unit



Closed position of lever tap outlet valve



Single-handed use

Flexible

Different tube lengths allow the user the choice of wearing the urine bag either on the thigh or on the lower leg; in three of the bags (Ref. 850353-000010, 850354-000000 and 850358-000000), the inlet tubing can even be cut to the desired length

Convenient

Clear graduation in ml to monitor the urine volume

Reliable

Flexible inlet tubing

Individual

Three-Chamber Leg Urine Collection Bags are available in both sterile and non-sterile versions

Caution

Non-sterile urine collectors should not be connected with sterile transurethral, suprapubic or nephrostomy catheters or any other sterile catheter. Non-sterile urine collectors can be connected only with non-sterile drainage devices (e.g. male external catheters).

Three-chamber Leg Urine Collection Bag

With lever tap outlet valve – Not made with natural rubber latex – Sterile – with 2 silicone backed comfort fabric straps per qty

CAPACITY	TUBE LENGTH	REF.	QTY
2000 ml	10 cm	850351-000000	8
	10 cm	850356-000000	10
	30 cm	850350-000000	8
	30 cm	850355-000000	10
	45 cm, variable	850353-000010*	10

With lever tap outlet valve – Not made with natural rubber latex – Non-sterile – with 2 silicone backed comfort fabric straps per qty

CAPACITY	TUBE LENGTH	REF.	QTY
500 ml	45 cm	850352-000000*	8
	45 cm	850357-000000	10
	45 cm	850359-000000***	20
	45 cm, variable	850354-000000*	8

With lever tap outlet valve – Not made with natural rubber latex – Sterile – with 2 pre-assembled standard button straps per bag

CAPACITY	TUBE LENGTH	REF.	QTY
750 ml	45 cm, variable	850358-000000*	50

* Without needle-free urine sampling port

** Individually packed in sealed transparent polybag