

# CERTIFICATE

APPROVED ALCOHOL SCREENING DEVICE

NHTSA/DOT/USCG



## -AL6000- TYPE TESTER

The model is a accurate and low cost alcoholtester.  
It is approved as an Alcohol Screening Device by: NHTSA - DOT - FDA510(k), and meets the US Coast Guard requirements.  
The digital readout is from 0 to 0,40% BAC, and special one-way mouthpieces are used to prevent false results.  
It has a built in counter to document the number of tests taken. Includes manual, batteries, carry case, 5 mouthpieces.

### SPECIFICATIONS

Alcohol sensor	: Semiconductor type alcohol sensor
Size	: 125 x 65 x 38 mm
Body material	: Plastic
Weight	: 117 g
Battery	: Two 1,5 V 'AA' size alkaline batteries
Detection range	: 0,00% - 0,40% BAC
Lifetime tests	: 3000 tests (~200 tests per batteries)
Recalibration need	: After placing new sensor the AL6000 auto-calibrates
Package	: Retail pack
Breath sampling time	: 4 seconds
Sensor accuracy	: (+/- 0,01 at 0,10% BAC)
Users	: Law enforcements
Warranty	: 1 year
Sensor	: 200 tests or 1 year

MONTH / YEAR / SIGNATURE



# CERTIFICATE

APPROVED ALCOHOL SCREENING DEVICE

NHTSA/DOT/USCG



## -AL7000- TYPE TESTER

The model is a accurate and low cost alcoholtester.  
It is approved as an Alcohol Screening Device by: NHTSA - DOT - FDA510(k), and meets the US Coast Guard requirements.  
The digital readout is from 0 to 0,40% BAC, and special one-way mouthpieces are used to prevent false results.  
It has a built in counter to document the number of tests taken. Includes manual, batteries, carry case, 5 mouthpieces.

### SPECIFICATIONS

Alcohol sensor	: Semiconductor type alcohol sensor
Size	: 105 x 50 x 18 mm
Body material	: Plastic
Weight	: 98 g
Battery	: Two 1,5 V 'AA' size alkaline batteries
Detection range	: 0,00% - 0,40% BAC
Lifetime tests	: 3000 tests (~200 tests per batteries)
Recalibration need	: After placing new sensor the AL7000 auto-calibrates
Package	: Retail pack
Breath sampling time	: 5 seconds
Sensor accuracy	: (+/- 0,01 at 0,10% BAC)
Users	: Law enforcements
Warranty	: 1 year
Sensor	: 200 tests or 1 year

MONTH / YEAR / SIGNATURE