

NOBAMED Paul Danz AG

NOBATEX[®] Individually packed in foil

REF 413360

Product Description, Intended use, Application

This white, woven bandage with a width of **6cm** and a **stretched length of 4m** is a lengthwise-elastic dressing retention bandage. Due to its weave, the bandage has a curled or creped structure. The elastic bandage is used to fix wound pads in which its padding effect is advantegous and is particularly suitable over joints and body curves. As a damp bandage, NOBATEX[®] is also used to shape and/or fix plaster splints. NOBATEX[®] is for single use.

Composition

Polyester, polyamide

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative

Requirements, Common Standards Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

NOBATEX[®] is manufactured according to DIN 61634.

Packaging

Primary packaging:	PP foil
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

Symbols used in labelling

Explanations at <u>www.nobamed.com</u>



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN 1041.

Storage and Transport

Dry and dustfree.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.