

TECHNICAL DATA SHEET

Vyconnector

801.00 D

Dated : 23-10-2024


1 Administrative information relating to the company

1.1	Name : Vygon	
1.2	Address : 5 rue Adeline - 95440 Ecouen, France	Fax : + 33 (0)1 34 29 19 34 E-mail : questions@vygon.com Website : www.vygon.com
1.3	Medical Device vigilance Representative : Laurent GUILLARDEAU	Tel : +33 (0)1 39 92 65 69 Fax : +33 (0)1 39 92 64 82 E-mail : quality@vygon.com


2 Information about the device or equipment

2.1					
2.2	Commercial name : Vyconnector				
2.3	GMDN code : 64467				
2.4	Medical Device Class : IIa Applicable regulation : 93/42/CEE in accordance with Appendix n° : II Notified body N° : 0459 Medical Device Manufacturer : Vygon				
2.5	Description of the device : Adaptor with a conical distal part and a female Luer-lock port. It allows the connection of a Luer syringe on the funnel of a feeding tube, vesical catheter etc. (for suctioning, lavage, drug injection...)				
					
2.6	Packaging / Containers				
	Code	Single unit packaging	Multi unit packaging	Minimum delivery quantity	Case
	801.00 D	1 (Supple blister)	100 (Carton box)	100 (Carton box)	2400 (Case box)

Single unit packaging



Box



Technical features :

Code	Connector
	Ø mm
801.00 D	7 - 10

2.7 Composition of the device and Accessories :

COMPONENTS	MARKER	MATERIALS
Connector		MBS



For components which are likely to come into contact with the patient and/or administered products, additional points :

Latex-free

DEHP-free

Does not contain any products of animal or biological origin

- 2.8 Indications :**
Connector to fit a Luer syringe to the funnel of a tube.

3 Accessories

4 Sterilization process

Sterile Medical Device : YES

Sterilization made for the device : Ethylene oxide

5 Conditions of conservation and storage

- 5.1 Normal conservation and storage conditions :**
Storage environment temperature: between 5 and 40°C. Store protected from moisture and sunlight.
- 5.2 Special precautions :**
Please refer to the insert in the packaging.
- 5.3 Duration of product validity :** 60 months

6 Security of use

Technical security :

Please refer to the insert in the packaging.

7 Instructions for use

- 7.1 Instructions :**
Please refer to the insert in the packaging.

7.2	Indications : Connector to fit a Luer syringe to the funnel of a tube.
7.3	Precautions : Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reaction or bacterial infections. Please refer to the insert in the packaging.
7.4	Contra indications : Please refer to the insert in the packaging.

8 Additional information relating to the product